

## Mouse anti-medaka Vtg Monoclonal antibody, CK-1H11

### General description

This product consists of affinity-purified mouse monoclonal antibodies against vitellogenin from Japanese medaka (*Oryzias latipes*).

Immunogen source: Lipovitellin purified from eggs of female Japanese medaka (*Oryzias latipes*).

Subclass/isotype: IgG<sub>1</sub>-kappa (and some IgG<sub>2b</sub>).

### Specificity

The monoclonal antibody CK-1H11 binds with high affinity to vitellogenin and lipovitellin from Japanese medaka (*Oryzias latipes*). The antibody also cross-reacts with vitellogenin from a variety of different species including Atlantic cod (*Gadus morhua*), rainbow trout (*Oncorhynchus mykiss*), carp (*Cyprinus carpio*) and tilapia (*Oreochromis niloticus*). The degree of cross-reactivity differs between species and with the methods employed.



### Applications

The monoclonal antibody CK-1H11 works well in both ELISA and western blot for a number of species. Since assay conditions vary, the optimum dilution should be determined for each particular application. *Note*: The antibody is dissolved in a buffer containing BSA, and is therefore not recommended to use for coating.

Normal dilution range:

ELISA: 1:100 - 1:1000

Western blot: 1:50 - 1:500



### ELISA

Coating: Purified vitellogenin (Vtg; 5 µg/ml) and lipovitellin (Vit; 5 µg/ml). Whole body homogenate (wbh; 1:1000) and plasma (pl; 1:1000) from female (F), 17β-estradiol treated (E2) or control (C) fish.

Primary antibody: CK-1H11

### Western blot

Samples: Purified vitellogenin (Vtg; 0.6 µg/well) and lipovitellin (Vit; 1 µg/well). Whole body homogenate (wbh; 1:20, 10 µl/well) from 17β-estradiol treated (E2) or control (C) fish.

Primary antibody: CK-1H11 diluted 1:100